


OPO



A P E



PP Automatic



The **OPO PP Automatic** is a synchronously pumped Optical Parametric Oscillator, designed for operation with modelocked Ti:Sapphire lasers. It is based on periodically poled crystals offering high conversion efficiency and wide tuning ranges with a fixed pump wavelength. Optical Parametric Oscillation is a nonlinear process converting a short wavelength pump beam into two tuneable beams (Signal + Idler) of longer wavelengths.

The exceptional tuning ranges of 1000-1600 nm (IR) and 505-750 nm (VIS-SHG) are covered with a single mirror set. An integrated spectrometer and the electronic tuning mechanism allow complete computer controlled operation of the OPO.

The **OPO PP Automatic** is supplied in different versions which are convertible and upgradeable.



Typical fs-OPO PP Automatic tuning curves, obtained with one optics set only at a fixed pump wavelength

Wide wavelength range

Computer controlled tuning and wavelength stabilization with integrated spectrometer

Femto- and picosecond mode

Low threshold

High conversion efficiency

Wavelength independent pulse duration

Low jitter - high stability

Ease of alignment

Ultrafast Pulse Diagnostics

Wavelength Conversion

Pulse Management


Acoustooptics

Your Partner in Ultrafast

OPO


SYNCHRONOUSLY PUMPED OPTICAL PARAMETRIC OSCILLATOR

The APE **OPO** covers a wide wavelength range and allows further wavelength extension using frequency mixing schemes.




CAVITY CONFIGURATIONS

The Linear IR - OPO version is a singly Signal resonant system with a 5-mirror standing wave cavity. It is used for highly efficient IR-generation and covers the 1,0...1,6 μm wavelength range and up to 3 μm with the non-resonant Idler branch (optional).



The Ring configuration is a singly Signal resonant IR/VIS **OPO** version employing an additional intracavity SHG module in an 8-mirror Ring cavity. This unit which is based on temperature tuned non-critically phasematched SHG is used for highly efficient visible output generation. It covers the 505 ... 750 nm wavelength range which fills the gap between the Ti:Sapphire Fundamental beam and its frequency doubled Signal (SHG). All Ring systems include the linear IR operation.



PULSE DURATION


All versions of the **OPO** system can be operated in both femto- and picosecond mode depending on pump pulse duration.

COMPONENTS

Control Electronics:



Optical Unit:



OPTIONS AND MODIFICATIONS

The APE **OPO** can be provided with several options and modifications:

- Autocorrelator
- Idler output (in addition to Signal output)
- Depleted pump output
- Customized wavelength ranges
- Adaptation to various pump laser cavity lengths and pump pulse durations
- Synchronization of different output beams (internal delay line) for use in frequency mixing devices and multicolor applications
- Frequency mixing devices (SumFrequencyMixing for extension to blue range, DifferenceFrequencyGeneration for Mid-IR generation)

SPECIFICATIONS

Tuning range (nm)	Pump	Signal (resonant)	Signal SHG ¹⁾	Idler (non-resonant)
Standard crystal PP830	830±10	1100 - 1600	550 - 750	1750 - 3300
Standard PP830+	830±10			
Extension crystal PP775	775±10	1000 - 1600	505 - 750	1750 - 3300
	Mira + Verdi 5W	Mira + Verdi 10W	Mira HP/Cham. Ultra II	
Signal output power @ 1300 nm	60 mW	180 mW	550 mW	
Ring SHG output power @ 600 nm	40 mW	150 mW	500 mW	
Threshold	500 mW ²⁾			
Pulse duration		typ. 200 fs @ 130 fs pump pulse typ. 1.6 ps @ 1.4 ps pump pulse		
Time-bandwidth product		typ. 0.6		
Polarization		horizontal for Signal and Idler, vertical for VIS-SHG		
Noise		<0,5 % RMS ³⁾		
Repetition rate		appr. 76 / 80 MHz, subject to pump laser repetition rate (others on request)		
Beam height		125 ... 140 mm		
Spectrometer	IR version	IR/VIS version	resolution	
	integrated into optics unit	external		
	700 ... 1640 nm	500 ... 1640 nm	0.2 nm	
Beam Parameters (typical values @ 1300nm)				
	Spatial mode	TEM ₀₀		
	M ²	1.2		
	Beam diameter	2 mm		
	Divergency	1 mrad		
Size (W*L*H in mm)	Optical unit	402 x 1139 x 200 (see p. 3)		
	Control electronics	267 x 312 x 180		

1) Ring-OPO with intracavity - SHG

2) low power systems; optimized for Mira + 5 W Verdi pump

3) RMS noise-measured in 10 Hz to 1 MHz bandwidth

Distributors

see APE website www.ape-berlin.com

APE GmbH Plauener Straße 163-165 Haus N / 13053 Berlin Germany

Phone +49.30.986.01130 Fax +49.30.986.011333 / Web www.ape-berlin.com Email ape@ape-berlin.de

APE follows a policy of continued product improvement. Therefore, specifications are subject to change without notice.