

Quantitative comparison study of flow cytometers using a novel ultra-stable calibration light source

Konrad von Volkmann, APE GmbH, Berlin

Oral Presentation, Emerging Technologies Session, 26th Annual Conference of the German Society for Cytometry.



Need for comparison between cytometers

- allow reproducible measurements
- at different times on the same instrument
- on different instruments (same site or at distinct sites)

T. Kalina *et al.*: EuroFlow standardization ..., Leukemia **26**, 1986–2010 (2012)

- System check
- Performance monitoring
- Automated data analysis
- others

Calibration standards

Standard ways today:

- Comparison is based on the measurements of latex microspheres
- Quality control, instrument characterization (Q & B), instrument setup ...
- Standardized „cells“, reproducible and stable

Limitations of today's standards

- Mainly relative calibration, absolute reference is demanding (variable gain, light collection, alignment)
- Variation in size due to the manufacturing process
intrinsic CV > 2%
- Variations from lot to lot
- Degradation process, aging: beads get lumpy

© 1992 Wiley-Liss, Inc.

Cytometry 13:822–830 (1992)

Noise, Sensitivity, and Resolution of Flow Cytometers¹

Harald B. Steen

Department of Biophysics, Institute for Cancer Research, The Norwegian Radium Hospital, 0310 Oslo, Norway

Received for publication January 8, 1992; accepted April 20, 1992

Cells passing through the laser focus = flashes of light

Detection of light by a PMT

- detection process is governed by Poissonian statistics
- intrinsic broadening of the detected mean value

$$CV = \frac{1}{\sqrt{N_{\text{photo electrons}}}} \quad \longleftrightarrow \quad N_{\text{photo electrons}} = \frac{1}{CV^2}$$

$$N_{\text{photo electrons}} = \text{Spe (statistical photoelectron estimate)}$$

© 1992 Wiley-Liss, Inc.


Cytometry 13:822-830 (1992)

Noise, Sensitivity, and Resolution of Flow Cytometers¹

Harald B. Steen

Department of Biophysics, Institute for Cancer Research, The Norwegian Radium Hospital, 0310 Oslo, Norway

Received for publication January 8, 1992; accepted April 20, 1992



www.quantiFlash.de

© 1992 Wiley-Liss, Inc.


Cytometry 13:822-830 (1992)

Noise, Sensitivity, and Resolution of Flow Cytometers¹

Harald B. Steen

Department of Biophysics, Institute for Cancer Research, The Norwegian Radium Hospital, 0310 Oslo, Norway

Received for publication January 8, 1992; accepted April 20, 1992



www.quantiFlash.de


Scale calibration in Spe (Statistical Photoelectron Unit)

quantiFlash® LED pulse -> Spe/Ch


Why Spe?

- Absolute reference
- Independent of gain
- Independent of filters used
- Independent of staining and biology

$$\text{Spe} = (100/rCV_{\text{pulse}})^2$$




Spe @ constant light intensity



Scale calibration in Spe

quantiFlash LED pulse -> **Spe/Ch**


$$\text{Spe} = (100/rCV_{\text{pulse}})^2$$



Scale calibration in ABC

measuring of beads with a calibrated
fluorescent dye value -> **ABC/Ch**

Quantum™ Simply Cellular® Microspheres




Calculation of Q

$$Q = \text{Spe} / \text{ABC} \quad (\text{or MESF})$$

Scale calibration in Spe

quantiFlash LED pulse -> Spe/Ch


$$\text{Spe} = (100/r\text{CV}_{\text{pulse}})^2$$



Sorter 1: FACSaria II SORP
lasers: 488, 561, ...


Sorter 2,3: FACSaria I
lasers: 488, ...

bandpass filter




In the PE Channel 585/20 and the PE-Cy7 Channel 780/60 the improved excitation efficiency leads to a 2-3 times higher detection efficiency.

quantiFlash® : schematic setup




Specifications of *quantiFlash*®

Pulse width	1-10 μ s, variable
Repetition rate	0.5-20 kHz, variable
Pulse shape	variable
Pulse amplitude	0 ... -96 dB
Pulse amplitude precision	CV < 0.1%
Fiber coupled	f-SMA termination
Power supply	rechargeable, USB powered




Graphical user interface for comfortable setup

- save/load parameter sets
- API for scripting (based on TCP/IP)
- easy to use software and firmware upgrades („one-click“)




Application examples using *quantiFlash*®

- Scale calibration (channels → Spe)
- Simulating multi-peak beads
- High-quality characterization measurement
- Optimized PMT voltages
- Check of detectors independent of laser alignment and sample preparation
- Check of linearity of detectors
- Study effects of pulse shape / pulse duration (check doublet discrimination)
- General alignment tool
- ...



Simulated 30-peak beads (3 dB intensity steps)
measured on FACSaria I, T. Kaiser / DRFZ

Calibrated data to Spe



1. Reference sample of single stained CD4-PE human PBMCs measured on two different sorters
2. Scale calibration to Spe using *quantiFlash*[®]


Sensitivity of sorter 1 is about 4× higher

Resolution is the relevant measure
depends on sensitivity **AND** background



Stain index


$$SI_1 = 22.3, SI_2 = 17.5$$



- Calibration is important for reproducible measurements and results
- *quantiFlash*®: a new tool for calibration
- Calibration of intensity scales leads to a meaningful basis for comparison



Calculation of Q
 $Q = \text{Spe} / \text{ABC} \quad (\text{or MESF})$



Thanks to...

Jan Popien, Sebastian Wolf, Ronny Schilling,
Frank Reuter & Katrin Keste



Toralf Kaiser, Kristen Feher, Jenny Kirsch &
Andreas Radbruch



Jim Wood,
Wake Forest University School of Medicine, USA



EUROPÄISCHE UNION
Europäischer Fonds für
regionale Entwicklung
Investition in Ihre Zukunft



**For further discussions visit
us at booth 16**

Q & B measurement (by Jim Wood)

Sample Name
Date

V570A
08-Dec-2015

Copyright, 2015; James C. S. Wood

Version 1.11 12/04/2015

CLEAR

PASTE

CALCULATE

Intrinsic CV

Spe per Channel

Background

Intensity Units

Spe

Background ERF

Goodness of Fit

Q

0.28%

6.16E-01

1.14E+02

7.03E+01

239.29

1.000

0.4770

